# **MONO-JET**

$\alpha$   $\beta$   $\gamma$   $\delta$

- MANUALE ISTRUZIONI
- OPERATOR'S HANDBOOK
- MANUEL D'UTILISATION
- GEBRAUCHSANWEISUNG
- MANUAL DE INSTRUCCIONES

**CATTANI S.p.A.**

# ***MONO-JET $\alpha$ $\beta$ $\gamma$ $\delta$***

## **INDICE**

## **Pagina**

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## **INDEX**

## **Page**

— FUNCTIONING GENERAL DATA .....	10
— LEGEND OF COMPONENTS .....	11
— SIGNALS AND WARNINGS .....	12
— GENERAL FEATURES .....	12
— INSTALLATION .....	12
— OPERATION AND USE .....	13
— GENERAL RECOMMENDATION AND BIOLOGICAL HAZARD .....	14
— MAINTENANCE AND CLEANING .....	14
— MAIN CLEANING AND MAINTENANCE OPERATIONS .....	15
— MAINTENANCE OPERATIONS MEANT FOR SERVICE ENGINEERS .....	15
— TRANSPORT AND STORAGE .....	16
— TRANSPORT OF SECOND-HAND APPLIANCES .....	16
— TABLE OF COMPONENTS .....	42

## Functioning general data

### Dental aspirator

<b>Model:</b>	<b>Mono-Jet <math>\alpha</math>, <math>\beta</math>, <math>\gamma</math>, and <math>\delta</math></b>
<b>Rated voltage:</b>	<b>230 V <math>\sim</math></b>
<b>Rated frequency:</b>	<b>50 Hz</b>
<b>Rated current:</b>	<b>3.1 A</b>
<b>Insulation class:</b>	<b>I</b>
<b>Type of apparatus:</b>	<b>B</b>
<b>Use:</b>	<b>continuous working</b>
<b>Level of protection against liquids:</b>	<b>common</b>
<b>Level of protection against direct or indirect contacts:</b>	<b>type B </b>

	<b>Alternating current</b>	<b>IEC 417-5032</b>
	<b>Protective earth (ground)</b>	<b>IEC 417-5019</b>
	<b>Apparatus type: B</b>	<b>IEC 878-02-02</b>
	<b>Off (power)</b>	<b>IEC 417-5008</b>
	<b>On (power)</b>	<b>IEC 417-5007</b>

**Mono-Jet  $\alpha$  and  $\beta$  230 V  $\sim$  50 Hz: authorization**

**N. J 0031**

**CEI 62-5**

**This appliance cannot work in the presence of an anaesthetic mixture inflammable with air, with oxygen or nitrogen protoxide.**

**Thermal motor protector.**

**Power output: 0.4 kW.**

**Maximum flow: 1250 l/min.**

**Maximum operating head for continuous service: 1300 mm H<sub>2</sub>O.**

**\*Sound pressure level 59 dB (A).**

**Other available tensions:**

**220 V  $\sim$  60 Hz 3.5 A**

**120 V  $\sim$  60 Hz 6.0 A**

**110 V  $\sim$  60 Hz 7.0 A**

**100 V  $\sim$  60 Hz 7.0 A**

**\*Sound pressure level tested (with canalised air) according to ISO 3746-1979 (E) regulation. Parameters: r or d = 1 - background noise  $\leq$  59 dB (A) - instruments: Brüel & Kjær Type 2232.**

**Manufactured by CATTANI S.P.A. - PARMA - ITALY**

## **LEGEND OF COMPONENTS** *(Table a page 42)*

- 1 - On/Off switch**
- 2 - Exhausted air silencer**
- 3A - Debris conveying pipe**
- 3B - Aspirated air conveying pipe to the motor**
- 3C - Exhausted air conveying pipe to silencer**
- 5 - Maximum level probes**
- 6 - Draining valve**
- 7 - Uni-Jet 75 aspiration unit**
- 8 - Maxi-Canister draining pump**
- 9 - Filter on canister cover**
- 10 - Debris tube-holder sleeve (on canister)**
- 11 - Maxi-Canister**
- 12 - Draining pump**
- 14 - Hydrocyclone**

## Signs and warnings

- **Biological danger, danger of infections from epidemic diseases.**


- **High temperature.**


- **General danger sign.**

- **Compulsory direction of flow and rotation.**

Signs cannot always fully express danger warnings, and the instructions deemed compulsory therefore it is necessary that the user reads the warnings and keeps them in due consideration. Failure to respect a danger sign or warning may harm operator or damage the equipment. Do not remove protections, do not tamper with machines or their operation.


ill. A

## General features (Mono-Jet inside view page 42)

**Mono-Jet is the aspirator designed to be coupled with the dental unit;** it is wieldy, noiseless and has been manufactured for continuous service.

Special versions:

- MONO-JET  $\alpha$  and  $\beta$ : differ only in the outside box.
- MONO-JET  $\gamma$ : draining is carried out during regular service.
- MONO-JET  $\gamma$  + HYDROCYCLONE: components as above, with Hydrocyclone.
- MONO-JET  $\delta$ : drainage by gravity, without draining pump.


ill. B

## Installation

The best position for the Mono-Jet is near the dental unit, precisely on the left of the chair, close to the patient's feet (ill. A).

**The aspirator must be installed in compliance with CEI regulations concerning electromedical appliances,** the same which have been complied with on designing and manufacturing our aspirators "Mono-Jet  $\alpha$ ,  $\beta$  and  $\gamma$ ", **CEI 62-5/2<sup>a</sup> 1991.** Before connecting the aspirator to the mains consult the ratings and **ascertain that the mains are suitable** for a proper and regular working of the appliance and that **it is protected from overcurrents and overloads as from CEI 64-4 regulation; the appliance must be installed in compliance with CEI 64-4** (relevant to appliances installed in medical rooms); **I.E.C. correspondent references are available on request for class I appliances.** The supply cable of the aspirator

ill. C

**is provided of protection cable for earthing;** this protection must not be removed for any reason and must fit perfectly in the supply plug. When the installation has been completed, the unit can be switched on by the main switch.

When the switch 1 is in ON position (ill. B), it will light up; lifting one of the terminals from its seat will start the aspiration. We recommend to use the brand-new, non contaminated appliance to train the surgery staff about the use and standard maintenance operations.

Opening of the door of the cabinet interrupts the circuit and the aspirator stops. On normal assembly the aspirated air is exhausted through the silencer 2 (page 42). Simply connect a pipe from the motor outlet to an outside tube 3C (ill. C), to convey the air outside: bacteria and most of the noise will be carried outside together with air.

## Operation and use

Through the tip and the lifted hose, liquids and aspirated air get to the manifold 4 (ill. H) and from there to the canister through pipe 3A. In the canister liquids are separated from air: the air through pipe 3B goes to the motor and is then exhausted; **on request it is available a certified bacteriological filter to filter the exhausted air.** Liquids, heavier than air, are collected on the bottom of the canister.

Draining is automatic: the valve 6 placed on the bottom of the canister is kept closed by depression when the aspirator is working; when all terminals are replaced on the tip support, the aspiration stops and liquids are drained from the canister (ill. D).

In case the canister should fill up during an operation, probes 5 (ill. D) will sense the maximum level at about  $\frac{3}{4}$  of the canister with the shortest probe; the electrical circuit will open and stop the aspiration automatically; the draining valve 6 will open automatically and the draining pump 8 (ill. E), will start: within few seconds the canister will be emptied and the aspiration can re-start automatically.

In the Mono-Jet type  $\gamma$ , fitted with draining pump, the aspiration is uninterrupted i.e. draining pump 8 (ill. F) starts without interrupting the aspiration.


**When using the Hydrocyclone (amalgam separator, 14 ill. F/2) consult manual "ISO Hydrocyclone and its applications".**

Blood, mucus and disinfectants can easily form foams that interfere with probes and interrupt suction.


The refore, use antifoam sanitizing tablets. If any running anomaly (block of the cooling system, breakdowns ecc.) causes an overheating of the motors i.e. temperature rises above 120 °C in the Uni Jet 75 (aspiration unit) and above 90 °C in the draining pump 8, in both cases thermal protectors with prefixed temperature value will open the circuits and reset them automatically as soon as the windings get back to normal temperature.

Obviously the causes of the interruption should be ascertained and rectified.


ill. D


ill. E

ill. F

ill. F/2

## General recommendations and biological hazard

- **Danger of electric shocks:**  
also 230 V ~ can lead to death.


- **Biological danger, danger of infections from epidemic diseases.**

Before any operation on used appliances, clean them several times with Puli-Jet or other adequate product.

**Cut off the appliance from the power mains and lock the general switch, if provided.**

**For any maintenance operation wear gloves, goggles, mask, and disposable overall. The aspirated liquid is infectious and can spread infective diseases.** In the Mono-Jet (ill. E and F, det. 8) the pump puts the aspirated fluid under pressure; for that reason, **the biological hazard increases if there is a pipe break, or a failure of the appliances under pressure.**

ill. G

## Maintenance and cleaning

- **Danger of electric shocks:**  
also 230 V ~ can lead to death.

- **Biological danger, danger of infections from epidemic diseases.**


Before operating on the appliance, read the chapter “**Signals and warnings**” and “**General notes and biological hazard**”


Because solid particles are aspirated together with liquids, restraining filters are fitted to assure a proper running. The Mini-Separator requires the filter 7 on the tip support (ill. H); the Mono-Jet fitted with Maxi-Canister has also a second filter 9 on the cover (ill. I). **Filters must be cleaned daily.** To reach the filter 7 the appliance must be switched on for some seconds with open terminals, aspirating air only so that hoses and manifold dry completely inside; then cut off the power, lift up the terminals from their seats and pull out the manifold 4 (ill. H): the handle of the filter will appear at the opening. To inspect the filter 9 lift up the pipe-holder sleeve 10 (ill. I). **The Maxi-Canister must be removed every week and cleaned with a jet of clear running water;** at the same time, the canister must be cleaned inside with a sponge as well as the lid and the overflow control probes. **Every evening, after cleaning the filters, it is very important that a solution of Puli-Jet and warm water should be aspirated; follow the instructions on the bottle label to prepare the solution.**

**This cleaning and sanitizing aspiration must be carried out by means of Pulse Cleaner (ill. L). With Pulse Cleaner the necessary turbulence for a complete cleansing is created.**

We strongly recommend avoiding the use of detergents, even with reduced foams, since turbulence and flow of aspirated air increase enormously the quantity of foams, causing interruptions, damage to the aspirator and in time also bad smells.

ill. H


ill. I

O-rings (seals) and the slides for the closure of terminals (ill. M) must be lubricated with spray silicon every 15 working days. It is important to replace all flexible tubes once a year (ill. P): in particular, the external tubes and terminals (ill. N) must be disinfected for each patient and replaced every three months for hygienic and operational reasons.

## Summary of main cleaning and maintenance operations

- After every surgical operation or very long operation: rinse the appliance by aspirating water, preferably warm. Replace and disinfect external tubes and terminals.
- After every working day: clean the filters; clean the equipment with warm water and Puli-Jet.
- Every 15 days: clean the canister, the draining valve and probes; lubricate O-rings and slides closing terminals with spray silicon.


ill. L

## Maintenance operations meant for service engineers


- Periodical checks: the service engineer shall check up all relays, test functioning of probes and draining pump, check siphons and waste outlets, all internal pipes, and plastic and rubber components subject to ageing.  
The appliance must be tested periodically from the current safety regulations point of view.

Short-circuit probes: make sure that suction stops and that the drainage pump starts operating. When using the Hydrocyclone, consult the manual **"ISO Hydrocyclone and its applications"**.

- The manufacturer is at your disposal to supply spare parts, literature, instructions and anything that might be useful. Our Concessionaires, Agents, authorised Retailers and Service Engineers are provided with split-up drawings, electrical diagrams, instructions and updating about maintenance and service.
- The appliance is guaranteed for a year from the date of sale provided that the guarantee card addressed to the manufacturer is returned with date of sale, retailer's stamp and customer's name.
- The guarantee and manufacturer's liability cease in case the appliances and/or the plants have been tampered with operations of any kind done by persons unfit and so unauthorised by the manufacturer.
- For any usage not contemplated or specified in this manual, please contact the manufacturer.
- On the web-site: [www.cattani.it](http://www.cattani.it) you can find our **updated** manuals. We suggest to consult them especially as for **safety** updating.

ill. M

ill. N

## ***Transport and storage***

- *Packed equipment can be transported and stored at a temperature range of -10 °C / +60 °C.*
- *Packages must be kept away from water and splashing and cannot tolerate humidity >70%.*
- *Packages with the same weight can be stored in piles of three only.*

## ***Transport of second-hand appliances***

- *Prior to packing cleanse and sanitise with PULI-JET or other suitable product (**see “maintenance and cleaning”**).*
- *Place unit into a polyethylene bag, seal and pack in 3-layer corrugated board.*

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